# SISSEL® BACK CHANGE

- Gebrauchsanweisung
- Instructions for use
- Mode d'emploi
- Gebruiksaanwijzing
- Instrucciones de uso
- Istruzioni per l'uso
- Bruksanvisning
- 사용 안내
- تعليمات الاستخدام

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

**A change for your back – our 2 in 1 back support and active sitting cushion**

- sitting pad with ergonomic design features a unique way to effectively guide the back between pressure and relief
- exercise and relax your back with ONE product:
  - for sitting: tension – active sitting improves core stability and works the deep muscles
  - placed behind your back: relief – perfectly supports the lumbar lordosis, enhances the natural sitting position and simultaneously incorporates back relaxation and gentle massage
- extreme comfort and usable by virtually everybody:
  - air-filled cushion with slightly rounded surface and straight bottom
  - approx. 40 cm diameter, load-bearing up to 100 kg
- degree of stability and firmness is freely adjustable by adding air (pump included)
- Polyester / Mesh cover for a proper air circulation (100 % polyester, Hand Wash)

**Application:**

**1. Correct sitting**

Large sitting surface (40 cm):

- for an improved pressure distribution and increased comfort while sitting
- ideal also for larger users

Adjustable seat height:

- by reducing some air it is possible to barely float above the sitting surface (needle pump included)

**2. Relaxation for your back:**

Back support and massage:

- Ideal supplement for dynamic sitting. Alternate between dynamic sitting, back relaxation respectively back massage.
- gentle, pleasant pressure on the back musculature
- supports the entire lower back (40 cm height): The air cushion contours to the individual shape of the spinal column and maintains the curvature of the lordosis.

**3. More than just a sitting cushion**

Wedge effect:

- The slight bulge enhances the physiological posture preventing the pelvis from sagging backwards and maintains a correct sitting posture.

No ridges at the edges:

- A smooth transition between cushion and chair surface without awkward edges.

Gentle rocking sensation at the rounded side:

- By reducing the air inside the cushion better meets the shape of the buttocks. The cushion reacts to even slight movements and transfers gentle vibrations into the deep musculature and this way strengthens core and back muscles.





**4. Attractive, 2-tone design**

- ideal to use at home, at work or when travelling
- this designer cushion complements your furniture
- air mesh cover ensures pleasant surface feel and air circulation

**5. Adjustable cushion firmness**

- The needle pump permits the addition or removal of air whenever desired. Just screw the tip of the needle onto the pump and carefully insert into the valve. Do this in a careful manner in order to avoid any damage of the SISSEL® Back Change outer membrane. Pressing your knee or lower arm onto the SISSEL® Back Change or sitting down on it is helpful because it creates a bulge around the valve (this avoids a puncturing of the product). Now insert the pump and add air. To deate simply the needle into the valve, hold slightly sideways until the desired amount of air has escaped.

**6. Use also for relaxation when resting on your back**

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**SISSEL® GmbH**